

Refrigerator drawers AluCool^{®} including dividers

High quality aluminium drawers including frame on smooth-running wheels with pull-out stop. Includes flexible dividing system with card pockets. Depending on the specific need, the dividing system allows the storage of drugs, blood and plasma with a space saving of up to 30 %. With a few simple steps, the flexible dividers are taken out and the drawer can be cleaned easily and hygienically.

The AluCool^{®} drawers can be retrofitted at any refrigerator - to save and optimise valuable space.

Tool-free assembly without damaging the refrigerator interior.

Division: 3 channels with 12 universal dividers. Width of the channels can be changed in 1cm steps.

For type:

FKUv 1610, 1613, 1660, 1663 / LKv 3910, 3913 / LKUexv 1610 / LKUv 1610, 1613 / LCv 4010 / GCv 4010, 4060 / MKv 3910, 3913 / MKUv 1610, 1613

Type	Description	PK	Cat. No.
B	For refrigerator internal dimensions 435 x 440 mm	1	9.698 981