

Molecular BioProducts™ RNase AWAY^{™} Surface Decontaminant

Eliminates unwanted DNA and RNase from glassware and plasticware without affecting subsequent DNA samples

Degrades DNA more quickly and effectively than autoclaving.

- Suitable for gel boxes, pipettes, benchtops, thermocyclers, or other apparatus
- Use before performing PCR or working with DNA for a DNA-free work area
- Simply apply the ready-to-use formula to the surface to decontaminate, then wipe dry or rinse clean

Use biocides safely. Always read the labelling and the product information before use.

Recommended application area: Surface disinfection

| Capacity ml | Supply form | Type of packaging | Bactericide | PK | Cat. No. |
|-------------|-----------------------|-------------------|-------------|----|-----------|
| 250 | Ready-to-use solution | Bottle | Yes | 1 | 6.225 974 |
| 475 | Ready-to-use solution | Spray flask | Yes | 1 | 6.227 799 |
| 1000 | Ready-to-use solution | Bottle | Yes | 1 | 6.206 320 |