

# Measuring pipette, serology, AR<sup>&reg;</sup>-glass, class B

Durable graduation with brown silver diffusion ink, wide opening. For serology. Graduation marks, with expanded tip opening, zero point at top, graduated to tip. With standard cotton plug from 5 ml.

Material: Soda-lime glass  
Colour graduation: Brown  
Accuracy class: B

Nominal capacity ml	Division ml	Tolerance ± ml	Length mm	With cotton plug end	PK	Cat. No.
1.0	0.01	0,01	360	No	1	6.225 600
2.0	0.02	0,015	360	No	1	6.225 601
5.0	0.01	0,045	360	Yes	1	6.056 581
10.0	0.1	0,075	360	Yes	1	6.510 390
25.0	0.1	0,15	450	Yes	1	6.510 389

Alternativen:

Graduated pipettes, serological, AR-GLAS<sup>®</sup>, total delivery / BRAND