

Bytac[®]: Bench protectors

Bytac[®]: Bench protectors, FEP Film with Vinyl or Aluminum Backing.

Inert to most chemicals - Non-stick; easy to clean - Use to resurface bench tops and line fume hoods. Contrasting work surface for viewing color changes during titrations. Protects surfaces against absorption of radioactive chemicals. Vinyl backing is recommended for all applications up to 93°C. Aluminium foil-backing is resistant to high temperatures up to 182°C for up to one hour.

Width mm	Length m	Material	Max. temp. resistance °C	Temperature resistance	PK	Cat. No.
635	4.6	FEP/Vinyl	93	93 °C	1	9.106 300
320	4.6	FEP/Vinyl	93	93 °C	1	6.267 458
50	33	FEP/Vinyl	93	93 °C	1	6.267 460
25	33	FEP/Vinyl	93	93 °C	1	6.267 461
635	4.6	FEP/Alu	182	182 °C	1	9.106 310
320	4.6	FEP/Alu	182	182 °C	1	6.225 629
50	33	FEP/Alu	182	182 °C	1	6.267 463