

Pipette controller accu-jet^{®} S, without power adapter

For precise pipetting of large and small volumes. The adjustable motor speed allows to switch between free flow and rapid blow-out of the liquid.

- One-hand operation
- For glass and plastic pipettes from 0.1 to 200 ml
- Precise volume setting
- Integrated rest position
- Valve system for exhaust liquid vapours
- Membrane filter and safety valve to protect against fluid penetration
- Robust, powerful motor for fast filling times (25 ml in 3 seconds)
- 8 hours of continuous pipetting with a 10 ml pipette without recharging
- Two-colour battery charge indicator
- Service-friendly and sustainable

Environmental benefits:

Sustainable packaging: Cardboard with approx 90 % recycled content
Renewable energy: Production using 100 % eco-friendly energy from certified hydro-electric power plants and energy generated on the premises as surplus from a cogeneration unit.

Scope of supply:

Pipette controller accu-jet^{®} S, nickel-metal hydride battery, wall support,
2 sterile spare membrane filters (0.2 μ m)
Power supply: 100 ... 240 V, 50/60 Hz

Min. capacity ml	Max. capacity ml	Colour	Type of auxiliary energy	PK	Cat. No.
0.1	200	Anthracite	Rechargeable battery	1	9.281 073
0.1	200	Berry	Rechargeable battery	1	9.281 070
0.1	200	Petrol	Rechargeable battery	1	9.281 071
0.1	200	Amethyst	Rechargeable battery	1	9.281 072

Zubehör:

Accessories for pipette controller accu-jet^{®} S / BRAND
Membrane filter for pipette controller accu-jet^{®} pro/S / BRAND