

Precision balances Explorer Plus EXP, with draftshield

The Explorer EXP series offers hassle-free data transfer with driverless HID (Human Interface Device). The time synchronization NTP function confirms data traceability to meet pharmaceutical regulations.

- 7 inch touchscreen
- Detachable terminal
- Motorised leveling system with digital leveling bubble
- Anti-pinch auto doors
- Touchless sensors
- Ionizers
- Full metal housing
- Removable draftshield
- Internal calibration systems
- 2 x USB-A, 1 x USB-B, 1 x USB-C, Ethernet, 1 x RJ11, 1 x RS232
- Optional: Bluetooth[®]/WiFi

Readability:	1 mg
Reproducibility:	0,7 mg
Linearity:	±0,6 mg
Weighing plate (WxD):	130 x 130 mm
With draft shield:	Yes
Dimensions (WxDxH):	229 x 408 x 372 mm
Weight:	7.8 kg
Power supply:	100 ... 240 VċAC, 50/60 Hz

Type	Max. load g	Verified calibration	Type of display	Width mm	Depth mm	Height mm	Type of auxiliary energy	Material hous
EXP223/AD	220		TFT	229	408	372	Mains connection	Metal
EXP223M/AD	220	Verified	TFT	229	408	372	Mains connection	Metal
EXP423/AD	420		TFT	229	408	372	Mains connection	Metal
EXP423M/AD	420	Verified	TFT	229	408	372	Mains connection	Metal
EXP623/AD	620		TFT	229	408	372	Mains connection	Metal
EXP623M/AD	620	Verified	TFT	229	408	372	Mains connection	Metal
EXP1203/AD	1200		TFT	229	408	372	Mains connection	Metal
EXP1203M/AD	1200	Verified	TFT	229	408	372	Mains connection	Metal

Alternativen:

- Precision balances Explorer Plus EXP / Ohaus
- Precision balances High Capacity Explorer Plus EXP / Ohaus
- Precision balances Explorer EXR, with draftshield / Ohaus
- Precision balances Explorer EXR / Ohaus