

Mixing cylinders, borosilicate glass 3.3, tall form, class A, blue graduated

All Blaubrand[®] cylinders are supplied with batch number and one batch certificate per packing unit. On request, they are also available with an individual certificate, USP individual certificate or DAkkS calibration certificate. Batch and individual certificates can also be downloaded at the supplier's website.

- Blaubrand[®]
- DIN EN ISO 4788
- DE-M marking
- PP stoppers and hexagonal foot
- Graduations and inscriptions in highly contrasting blue enamel
- Calibrated to contain (IN)

Height: 395 mm
Ground size: NS34/35
Type of closure: Stopper

Nominal capacity ml	Division ml	Shape	Accuracy class	Colour graduation	Material closure	Material	PK	Cat. No.
500	5.0	Tall	A	Blue	PE	Boro 3.3	1	9.275 254

Alternativen:

Mixing cylinders, borosilicate glass 3.3, stoppered, class B, amber graduated / BRAND
Mixing cylinders, borosilicate glass 3.3, tall form, class A, blue graduated / BRAND