

# Refrigerator FKU 1800

The refrigerator can be integrated underneath a worktop to save space.

- Efficient refrigeration system with static cooling
- Digital temperature display
- LED interior lighting
- Automatic defrosting
- Lockable door
- Door hinge right, exchangeable
- Environmentally friendly, natural refrigerant R 600a

Temperature: 2 ... 10 &deg;C  
Ext. dim. (WxDxH): 600 x 600 x 850 mm  
Int. Dim. (WxDxH): 513 x 441 x 702 mm  
No. of chambers: 1  
No. of shelves: 3  
No. of doors: 1  
With glass door: No  
Door hinge: Changeable  
Weight: 33 kg  
Energy consumption (24h): 0.559 kWh  
Power supply: 220 ... 240 V, 50 Hz

Type	Nominal capacity l	Min. temperature °C	Max. temperature °C	Ext. width mm	Ext. depth mm	Ext. height mm	Int. width mm	D m
FKU 1800	160	2	10	600	600	850	513	44

Zubehör:  
Refrigerator drawers AluCool® including dividers / H+H System