

Filter tips, PP

The pipette tips are produced in clean room and are compatible with most popular brands of pipettors.

- Sterile
- DNase-, RNase- and pyrogen free
- Graduated

Colour: Transparent
Tip form: Standard
With graduation: Yes
Package: 10 x Rack of 96
Sterile: Yes

Capacity	Min. capacity µl	Max. capacity µl	Material	RNase/DNase-free	Pyrogen free	PK	Cat. No.
0.5 ... 10 µl	0.5	10.0	PP	Yes	Yes	960	4.682 025
1.0 ... 20 µl	1.0	20.0	PP	Yes	Yes	960	4.682 046
10 ... 100 µl	10.0	100.0	PP	Yes	Yes	960	4.682 047
10 ... 200 µl	10.0	200.0	PP	Yes	Yes	960	4.682 048

Alternativen:
LLG-Filter tips 2.0 Low Retention, PP / LLG Labware