

Disposable Syringes Omnifix[®] Solo, 3-piece, Luer Lock conus

- Material: Polypropylene, Plunger backstop: Polyisoprene
- Plunger backstop with smooth double sealing ring for slow aspiration and injection of small sample quantities
- Highly transparent barrel with black graduations (in ml) for ideal readability, permanent marking
- Safe plunger backstop for easy aspiration to the maximum volume
- Inside siliconized
- Luer Lock conus, centric
- Latex- and PVC-free
- Sterile, single packed
- Manufactured according to EN ISO 7886-1

Nozzle type: Luer lock
Sterile: Yes
Package: Single packed

Nominal capacity ml	Material barrel	Material plunger	Colour graduation	Division ml	Standard	PK	Cat. No.
3	PP	PP	Black	0.1	EN ISO 7886-1	100	6.238 514
5	PP	PP	Black	0.2	EN ISO 7886-1	100	6.081 232
10	PP	PP	Black	0.5	EN ISO 7886-1	100	6.085 753
20	PP	PP	Black	1	EN ISO 7886-1	100	6.083 393
30	PP	PP	Black	1	EN ISO 7886-1	100	6.300 903
50	PP	PP	Black	1	EN ISO 7886-1	100	4.665 943

Alternativen:
Fine Dosage Syringes Omnifix[®]-F, 3-piece / B. Braun Deutschland
Disposable Syringes Omnifix[®] Solo, 3-piece, Luer slip / B. Braun Deutschland