

QuikSip®; BT-aspirator

The bottletop aspirator is designed for safe and fast aspiration of common laboratory liquids, e.g. biological solutions, nutrient media, polar solvents, aqueous solutions.

- Safe removal of supernatants (up to max. 25 ml per plunger stroke)
- Works without a vacuum pump
- Vacuum control using cell-culture® unit
- Works as single channel or 8-channel unit (optional)
- Suitable for use with disposable pipette tips, micropipettes and glass pasteur pipettes
- Adapter and suction tube for the cell-culture® unit are autoclavable (dispensing cartridge and pump unit are not autoclavable)

Scope of supply: 1 QuikSip®; BT aspirator, spare dispensing cartridge, 2 PP adapters (GL 45/32 and GL 45/S40)

Please order cell-culture®-unit separately.

Description	PK	Cat. No.
QuikSip®; BT-Aspirator	1	4.684 523

Zubehör:
Accessories for seripettor® pro / QuikSip™ BT-aspirator / BRAND

Alternativen:
Safety suction system EcoVac / schuett-biotec