

pH combination electrodes InPro 3250

Low maintenance pH electrodes with pre-pressurised liquid electrolytes for use in biotechnology and pharmaceutical industry. With integrated temperature sensor for temperature compensation.

- Sensor and plug completely autoclavable and sterilisable
- Approved for Ex-protected environments
- Certified in accordance with the DGRL
- With silver ion barrier for the measurement of sulfide-containing samples
- Pressure stability up to 4 bar
- "i" version with integrated sensor management ISM
- Reference system: Argenthal®

Please order connection cable separately.

pH range:	0 ... 12
Temperature:	0 ... 100 °C
Diaphragm:	Ceramic
Shape membrane:	round
Shaft Ø:	12 mm

Type	Max. pH value	Max. temperature °C	Temperature sensor	Connection plug	Shaft length mm	Material shaft
InPro 3253/150/PT100	12	100	Pt 100	VP plug head	150	Glass
InPro 3253i/SG/120	12	100	NTC	K8S	120	Glass
InPro 3253i/SG/225	12	100	NTC	K8S	225	Glass
InPro 3253/SG/120/PT1000	12	100	Pt 1000	VP plug head	120	Glass