

Multi-layer cell culture systems CellSafe™ CellFac®reg, 5-layer

- ISO 13485 and ISO 9001 certified
- Produced in class 100000 clean rooms
- Made with USP Class VI standards medical-grade raw materials
- GMP compliant
- Marked with batch number
- Triple-bagged packaging
- Sterilised by irradiation, SAL 10⁻⁶
- DNase/RNase-free, pyrogen-free, non-cytotoxic

Coating: TC
Min. working volume: 650 ml
Max. working volume: 1000 ml

Working volume	Growth surface cm ²	Sterile	RNase/DNase-free	Pyrogen free	PK	Cat. No.
650 ... 1,000 ml	3216	Yes	Yes	Yes	4	4.714 449