

# Multi-layer cell culture systems CellSafe™ CellFac®<sup>reg</sup>, 10-layer

- ISO 13485 and ISO 9001 certified
- Produced in class 100000 clean rooms
- Made with USP Class VI standards medical-grade raw materials
- GMP compliant
- Marked with batch number
- Triple-bagged packaging
- Sterilised by irradiation, SAL 10<sup>-6</sup>
- DNase/RNase-free, pyrogen-free, non-cytotoxic

Coating: TC  
Min. working volume: 1300 ml  
Max. working volume: 2000 ml

Working volume	Growth surface cm <sup>2</sup>	Sterile	RNase/DNase-free	Pyrogen free	PK	Cat. No.
1,300 ... 2,000 ml	6416	Yes	Yes	Yes	2	4.714 450

Alternativen:  
Multi-layer cell culture systems CellSafe™ CellFac®, 5-layer / Guangzhou JET Bio-Filtration