

# Cleanroom Gloves, BioClean ULTIMATE™, polychloroprene, sterile

- Low particulate count
- Powder-free
- Latex-free
- AQL 0.65
- Textured
- Hand specific
- Excellent ESD properties
- Easy double-donning
- Chemical resistant
- Compatibility ISO Class 4
- Colour: natural
- Carton with 20 bags with 10 pairs

According to EN ISO 374-1:2016/Type C, EN ISO 374-5:2016, EN 421, Suitable for food handling.

Packaging: 1 pair per plastic bag, 1 plastic bag per particle-free EasyTear™ bag, 10 pairs per sealed outer bag, 20 outer bags per inner carton bag

Colour: Natural  
Sterile: Yes  
Material: Neoprene  
Powder free: Yes  
AQL: 0.65

Glove size	Length mm	Suitable for cleanroom	Food safe	Ambidextrous	PK	Cat. No.
6	300	Yes	Yes	No	200	4.661 444
6.5	300	Yes	Yes	No	200	4.661 445
7	300	Yes	Yes	No	200	4.661 446
7.5	300	Yes	Yes	No	200	4.661 447
8	300	Yes	Yes	No	200	4.661 448
8.5	300	Yes	Yes	No	200	4.661 449
9	300	Yes	Yes	No	200	4.661 450