

Somatic cell counter ADAM™ SCC 2

ADAM™ SCC 2 is a somatic cell counter integrated with fluorescence optic and image analysis software for milk quality control.

- Built-in tablet
- Less than 13 seconds for a test (4 channel)
- PI staining method
- No need for calibration
- WiFi (USB Dongle)

Dimensions (WxDxH): 276 x 227 x 270 mm
Weight: 7 kg
Power supply: 12 VċDC
Warranty: 1 year

Type	Max. volume sample μl	Width mm	Depth mm	Height mm	Type of auxiliary energy	PK	Cat. No.
ADAM™ SCC 2	13	276	227	270	Mains connection	1	4.684 382

Zubehör:
SOMA Chip Kits / NanoEn Tek Inc.