

# Cleanroom Gloves, BioClean NERVA, Nitrile

- Elbow length protection
- Low particulate count
- Powder-free
- Latex-free
- AQL 0.65
- Textured
- Ambidextrous
- Excellent ESD properties
- Easy double-donning
- Chemical resistant
- Compatibility ISO Class 4
- Colour: white
- Carton with 10 bags with 100 pieces

According to EN ISO 374-1:2016/Type C, EN ISO 374-5:2016, EN 421, Suitable for food handling.

Packaging: 100 pieces per sealed plastic inner bag, 1 inner bag per sealed outer bag, 10 outer bags per inner carton bag

Colour: White  
Material: Nitrile  
Powder free: Yes  
AQL: 0.65

Glove size	Length mm	Suitable for cleanroom	Food safe	Ambidextrous	PK	Cat. No.
S (6.5 - 7)	400	Yes	Yes	Yes	1000	4.661 387
M (7.5 - 8)	400	Yes	Yes	Yes	1000	4.661 388
L (8.5 - 9)	400	Yes	Yes	Yes	1000	4.661 389
XL (9.5 - 10)	400	Yes	Yes	Yes	1000	4.661 390