

Cleanroom Gloves, BioClean ADVANCE®; latex, sterile

- Low particulate count
- Powder-free
- AQL 0.65
- Textured
- Hand specific
- Easy double-donning
- Chemical resistant
- Compatibility ISO Class 4
- Colour: natural
- Carton with 20 bags with 10 pairs

According to EN ISO 374-1:2016/Type C, EN ISO 374-5:2016, EN 421, Suitable for food handling.

Packaging: 1 pair per plastic bag, 1 plastic bag per particle-free EasyTear® bag, 10 pairs per sealed outer bag, 20 outer bags per inner carton bag

Colour: Natural
Sterile: Yes
Material: Latex
Powder free: Yes
AQL: 0.65

Glove size	Length mm	Suitable for cleanroom	Food safe	Ambidextrous	PK	Cat. No.
6	295	Yes	Yes	No	200	4.661 396
6.5	295	Yes	Yes	No	200	4.661 397
7	295	Yes	Yes	No	200	4.661 398
7.5	295	Yes	Yes	No	200	4.661 399
8	295	Yes	Yes	No	200	4.661 400
8.5	295	Yes	Yes	No	200	4.661 401
9	295	Yes	Yes	No	200	4.661 402
10	295	Yes	Yes	No	200	4.661 403