

Cleanroom Gloves, BioClean SYNERGY™, Nitrile

- Super-sensitive and accelerator-free
- Ambidextrous
- Powder-free
- Latex-free
- AQL 1.5
- Textured
- Excellent ESD properties
- Easy double-donning
- Chemical resistant
- Compatibility ISO Class 4
- Colour: green
- Carton with 10 bags with 100 pieces

According to EN ISO 374-1:2016/Type C, EN ISO 374-5:2016, EN 421, Suitable for food handling.

Packaging: 100 pieces per sealed plastic inner bag, 1 inner bag per sealed outer bag, 10 outer bags per inner carton bag

Colour: Green
Sterile: No
Material: Nitrile
Powder free: Yes
AQL: 1.5

Glove size	Length mm	Suitable for cleanroom	Food safe	Ambidextrous	PK	Cat. No.
XS (5.5 - 6)	300	Yes	Yes	Yes	1000	4.661 391
S (6.5 - 7)	300	Yes	Yes	Yes	1000	4.661 392
M (7.5 - 8)	300	Yes	Yes	Yes	1000	4.661 393
L (8.5 - 9)	300	Yes	Yes	Yes	1000	4.661 394
XL (9.5 - 10)	300	Yes	Yes	Yes	1000	4.661 395