

Indicator tape for steam, hot air and gas sterilisation

Crepe paper, adhesive. With indicator. Colour changes according to DIN EN ISO 11140-1.

Rolls of 50 m.

Length: 50 m
Content: 1

Colour change	Width mm	No. of packages	PK	Cat. No.
Green/black	19	1	1	9.140 784
White/black	25	1	1	9.140 789
Green/black	19	1	1	9.140 790
White/black	25	1	1	9.140 791

Alternativen:

Indicator Tape, Comply™ / 3M Deutschland