

# Cleanroom Gloves, BioClean FUSION®;, polychloroprene, sterile

- Ambidextrous
- Powder-free
- Latex-free
- Finger-textured
- Excellent ESD properties
- Easy double-donning
- Chemical resistant
- Compatibility ISO Class 4
- Carton with 20 bags with 10 pairs

According to EN ISO 374-1:2016/Type C, EN ISO 374-5:2016, EN 421, Suitable for food handling.

Packaging: 1 pair per plastic bag, 1 plastic bag per particle-free EasyTear®; bag, 10 pairs per sealed outer bag, 20 outer bags per inner carton bag

Colour: Green  
Sterile: Yes  
Material: Neoprene  
Powder free: Yes  
AQL: 0.65

Glove size	Length mm	Suitable for cleanroom	Food safe	Ambidextrous	PK	Cat. No.
5.5	300	Yes	Yes	Yes	200	4.661 451
6	300	Yes	Yes	Yes	200	4.661 452
7	300	Yes	Yes	Yes	200	4.661 453
8	300	Yes	Yes	Yes	200	4.661 454
9	300	Yes	Yes	Yes	200	4.661 455
10	300	Yes	Yes	Yes	200	4.661 456