

pH combination electrodes InPro 3250

pH electrodes with pre-pressurised liquid electrolytes for use in chemical processes. With integrated temperature sensor for temperature compensation.

- Low maintenance
- Approved for Ex-protected environments
- Certified in accordance with the DGRL
- With silver ion barrier for the measurement of sulfide-containing samples
- Pressure stability up to 4 bar
- Ag/AgCl reference system.

Please order connection cable separately.

pH range: 14
Temperature: 0 ... 100 °C
Shape membrane: round
Shaft Ø: 12 mm
Shaft length: 120 mm

Type	Max. pH value	Max. temperature °C	Temperature sensor	Connection plug	Material shaft	PK	Cat. No.
InPro 3250/120/PT1000	14	100	Pt 1000	VP plug head	Glass	1	6.280 818
InPro 3250i/SG/120	14	100	Pt 100	K8S	Glass	1	6.242 095