

# Strip plates 12 x F8 BRANDplates<sup>&reg;</sup>, immunoGrade&trade;, with grid

High binding capacity (highbinding), ideal for IgG and molecules with hydrophobic and hydrophilic regions.

12 strips of 8 wells each. Individual strips are removable from the plate and can be reinserted. The strips can be divided into 8 individual vessels and can be inserted individually into the frame plate.

Material: PS  
Colour: Transparent  
Bottom shape: Flat

No. of wells	Nominal capacity $\mu$ l	PK	Cat. No.
96	350	100	4.007 913