

Wireless module IDS WLM-S

Wireless measurement of pH/Redox, conductivity, dissolved oxygen and turbidity.

- Bluetooth^{®} LE
- Universal plug head
- Opening for electrolyte solution refilling (only pH electrodes with liquid electrolyte)
- Signal LED (RGB)
- Range 10 m
- IP code IP 66
- Power supply: Re-chargeable LiPo battery 230 mAh

Suitable for IDS portable meters MultiLine^{®}: Multi 3630 IDS/3620
IDS/3510 IDS, as well as for IDS laboratory meters inoLab^{®}: Multi 9630
IDS/9620 IDS/9310 IDS.

Type	PK	Cat. No.
IDS WLM-S	1	6.283 617

Zubehör:
USB charger for wireless module IDS WLM-S / Xylem Analytics Germany (WTW)