

Magnetic stirrer C-MAG HS 7

Connection acc. to DIN 12878 for attaching a contact thermometer, e.g. ETS-D5; enables precise temperature control.

Max. stirring capacity: 10 l
Speed: 100 ... 1,500 1/min
Heating temperature: 50 ... 500 °C
Max. heating rate: 5.0 °C/min
Working area (WxD): 180 x 180 mm
Material stirring surface: Glass ceramic
Colour housing: Blue
Safety temperature: 550 °C
Permitted ambient temperature: 5 ... 40 °C
Dimensions (WxDxH): 220 x 330 x 105 mm
Weight: 5 kg
Heating capacity: 1000 W
Rated capacity: 270 W
Power supply: 230 V, 50/60 Hz
IP code: IP 21

Type	No. of stirring positions	Min. speed 1/min	Max. speed 1/min	Min. heating temperature °C	Max. heating temperature °C
C-MAG HS 7	1	100	1500	50	500
C-MAG HS 7	1	100	1500	50	500
C-MAG HS 7	1	100	1500	50	500

Zubehör:

Magnetic stirring bars Ikaflon®, cylindrical, PTFE / IKA-Werke
Magnetic stirring bar set RS 1 / RS 2 / IKA-Werke

Alternativen:

Magnetic stirrers C-MAG MS 4 / C-MAG MS 7 / C-MAG MS 10 / IKA-Werke
Magnetic stirrer C-MAG HS 10 / IKA-Werke
Magnetic stirrer C-MAG HS 4 / IKA-Werke