

Tabletop autoclaves CertoClav Classic

CertoClav Classic is suitable for sterilisation of unwrapped solids without hollow bodies, liquids, bottles, instruments, synthetic materials, culture media preparation and waste sterilisation.

A precision thermometer, a control valve and a precision manometer with drag indicator ensure exact processes and perfect results for all applications. The model is available in two versions, each with two temperature levels. The respective temperature level is set via the control valve in the lid. The valves cannot be interchanged as the thermostat of the appliance is matched and adjusted to the respective valve. The water supply is a mixture of tap water and demineralised water.

Our autoclaves are protected by a multiple safety system. It consists of a pressure relief valve, a pressure valve, a safety lid lock to prevent opening under pressure and a BI metal fuse to prevent the appliance from overheating.

- Touchscreen control
- Robust construction, low-maintenance
- Transportable and location-independent
- Thermostat control
- Integrated timer
- Venting in flow process
- Heating integrated in the base

Scope of supply: Autoclave, instrument tray, power cable

Max. pressure:	2.7 bar
Material interior:	Stainless steel
Dim. Int.(Ø x H):	260 x 350 mm
Ext. dim. (WxDxH):	445 x 300 x 605 mm
Weight:	8.8 kg
Heating capacity:	1300 W
Power supply:	230 V, 50/60 Hz

Type	Nominal capacity l	Temperature levels	Int. Ø mm	Height int. mm	Ext. width mm	Ext. depth mm	Ext. height mm	Type of auxiliary energy
Classic	18	125/140 °C	260	350	445	300	605	Mains connection
Classic	18	115/121 °C	260	350	445	300	605	Mains connection

Zubehör:
Accessories for autoclaves, CertoClav / CertoClav Sterilizer

Alternativen:
Tabletop autoclaves CertoClav MultiControl / CertoClav Sterilizer