

Analytical balances Explorer Plus EXP

The Explorer EXP series offers hassle-free data transfer with driverless HID (Human Interface Device). The time synchronization NTP function confirms data traceability to meet pharmaceutical regulations.

- 7 inch touchscreen
- Detachable terminal
- Motorised leveling system with digital leveling bubble
- Anti-pinch auto doors
- Touchless sensors
- Ionizers
- Full metal housing
- Removable draftshield
- Internal calibration systems
- 2 x USB-A, 1 x USB-B, 1 x USB-C, Ethernet, 1 x RJ11, 1 x RS232
- Optional: Bluetooth[®]/WiFi

Readability:	0,1 mg
Reproducibility:	0,05 mg
Linearity:	±0,06 mg
Weighing plate (WxD):	90 x 90 mm
Dimensions (WxDxH):	229 x 408 x 372 mm
Weight:	7.8 kg
Power supply:	100 ... 240 V·AC, 50/60 Hz

Type	Max. load g	Verified calibration	Type of display	With draft shield	Width mm	Depth mm	Height mm	Type of auxiliary energy
EXP124/AD	120		TFT	Yes	229	408	372	Mains connection
EXP124M/AD	120	Verified	TFT	Yes	229	408	372	Mains connection
EXP224/AD	220		TFT	Yes	229	408	372	Mains connection
EXP224M/AD	220	Verified	TFT	Yes	229	408	372	Mains connection
EXP324/AD	320		TFT	Yes	229	408	372	Mains connection
EXP324M/AD	320	Verified	TFT	Yes	229	408	372	Mains connection