

Tabletop autoclaves CertoClav EL

The vertical laboratory autoclave is suitable for sterilisation of unpackaged solids without hollow bodies, liquids, bottles, instruments, synthetic materials, food preservation, culture media preparation and waste. A precision thermometer, a control valve, and a precision manometer with drag indicator provide precise processes and accurate results for all applications. The autoclave is available in two versions, each with two temperature levels. These temperature settings are controlled via the control valve on the lid. These valves are not interchangeable, because the thermostat of the device is harmonized with and calibrated to the respective valve.

- Robust construction
- Thermostat controlled
- High heating power
- Venting by flow procedure
- Heating integrated in the base

Scope of supply: Autoclave, instrument tray, power cable
Max. pressure: 2.7 bar
Material interior: Aluminium
Heating capacity: 1900 W
Power supply: 230 V, 50/60 Hz

Type	Nominal capacity l	Temperature levels	Int. Ø mm	Height int. mm	Dim. Int. (Ø x H)	Ext. width mm	Ext. depth mm	Ext. height mm
EL 12	12	125/140 °C	240	260	240 x 260 mm	395	275	505
EL 12	12	115/121 °C	240	260	240 x 260 mm	395	275	505
EL 18	18	125/140 °C	240	390	240 x 390 mm	395	275	665
EL 18	18	115/121 °C	240	390	240 x 390 mm	395	275	665

Zubehör:
Accessories for autoclaves, CertoClav / CertoClav Sterilizer

Alternativen:
Tabletop autoclaves CertoClav MultiControl / CertoClav Sterilizer