

Flask Inlay Series H 135 for IKA Magnetic Stirrers

For holding round flasks when used with magnetic stirrers. Made of anodised aluminium.
Can be used up to max. 180 °C.

Fits on magnetic stirrers: RCT basic, RET basic, RH basic, RH basic white, RH digital, RH digital white, RET control-visc, RET control-visc white, IKA plate (RCT digital), C-MAG HS 7 control (only with adapter).

Please order the appropriate inlay for the designated flask size.

Type	For	For vessels ml	PK	Cat. No.
H 135.201	H 135.20/H 135.21	10	1	9.720 823
H 135.202	H 135.20/H 135.21	25	1	9.720 824
H 135.203	H 135.20/H 135.21	50	1	9.720 825
H 135.301	H 135.30/H 135.31	100	1	9.720 828
H 135.302	H 135.30/H 135.31	250	1	9.720 829
H 135.401	H 135.40/H 135.41	500	1	9.720 832
H 135.501	H 135.50/H 135.51	1000	1	9.720 845