

Temperature chamber ETD 300 for rotational viscometers VisQC

- Temperature control from +25 to +300 °C
- Full control via the ViscoQC display
- For measurement of bitumen, hot melts, polymers, waxes
- Complies with ASTM D4402, D3236, AASHTO T316, EN 13302 and more

Type	For	PK	Cat. No.
ETD 300	ViscoQC 100, ViscoQC 300	1	4.703 002

Zubehör:
Cooling adapter for temperature chamber ETD 300 / Anton Paar Germany
Spindles for temperature chamber ETD 300 / Anton Paar Germany