

Microcentrifuge Frontier™; FC5513 IVD with angle rotor R02

Microcentrifuge with angle rotor 24 x 1.5/2.0 ml BIOSEALS

Speed: 200 ... 13,500 1/min
With cooling: No
Noise level: 60 db(A)
Dimensions (WxDxH): 277 x 356 x 236 mm
Weight: 13.1 kg
Permitted ambient temperature: 5 ... 35 °C
Max. relative humidity rate: 80 %
Rated capacity: 170 W
Power supply: 230 VA, 50/60 Hz

Type	Max. rcf x-g	Min. speed 1/min	Max. speed 1/min	Width mm	Depth mm	Height mm	Min. permitted ambient temperature °C
FC5513+R02 IVD	17317	200	13500	277	356	236	5

Zubehör:
Fixed angle rotors for micro centrifuges Frontier™ 5000 Micro / Ohaus

Alternativen:
Microcentrifuge Frontier™ FC5513 IVD / FC5513R IVD, without rotor / Ohaus
Microcentrifuge Frontier™ FC5513 IVD with angle rotor R06 / Ohaus