

# Volumetric flasks, boro 3.3, class A, blue graduations, incl. ISO individual certificate

All Blaubrand<sup>®</sup> volumetric flasks are delivered with batch number and one batch certificate included per pack. They are also available with DAkkS calibration certificate upon request.

- DIN EN ISO 1042
- DE-M marking
- Ring marks and inscriptions in highly contrasting blue enamel
- Calibrated to contain (IN)
- With PP- / PE-stoppers

Standard: DIN EN ISO 1042

Nominal capacity ml	Tolerance ml	Ground size	Accuracy class	Material	Material closure	Type of closure	PK	Cat. No.
5	&plusmn;0.04	NS10/19	A	Boro 3.3	PP	Stopper	1	7.659 411
10	&plusmn;0.025	NS7/16	A	Boro 3.3	PP	Stopper	1	4.684 479
10	&plusmn;0.04	NS10/19	A	Boro 3.3	PP	Stopper	1	7.659 412
20	&plusmn;0.04	NS10/19	A	Boro 3.3	PP	Stopper	1	7.085 453
20	&plusmn;0.06	NS12/21	A	Boro 3.3	PP	Stopper	1	4.684 480
25	&plusmn;0.04	NS10/19	A	Boro 3.3	PP	Stopper	1	7.659 413
25	&plusmn;0.06	NS12/21	A	Boro 3.3	PP	Stopper	1	4.684 481
50	&plusmn;0.06	NS12/21	A	Boro 3.3	PP	Stopper	1	6.253 094
50	&plusmn;0.10	NS14/23	A	Boro 3.3	PP	Stopper	1	6.257 937
100	&plusmn;0.10	NS12/21	A	Boro 3.3	PP	Stopper	1	4.684 482
100	&plusmn;0.10	NS14/23	A	Boro 3.3	PP	Stopper	1	6.256 708
150	&plusmn;0.15	NS14/23	A	Boro 3.3	PP	Stopper	1	4.684 483
200	&plusmn;0.15	NS14/23	A	Boro 3.3	PP	Stopper	1	7.400 043
250	&plusmn;0.15	NS14/23	A	Boro 3.3	PP	Stopper	1	6.264 670
500	&plusmn;0.25	NS19/26	A	Boro 3.3	PP	Stopper	1	6.264 671
1000	&plusmn;0.4	NS24/29	A	Boro 3.3	PP	Stopper	1	6.264 672
1000	&plusmn;0.6	NS29/32	A	Boro 3.3	PP	Stopper	1	4.684 484
2000	&plusmn;0.6	NS29/32	A	Boro 3.3	PP	Stopper	1	6.264 673
5000	&plusmn;1.2	NS34/35	A	Boro 3.3	PE	Stopper	1	4.684 485
10000	&plusmn;2.0	NS45/40	A	Boro 3.3	PE	Stopper	1	4.684 486

## Alternativen:

- Volumetric flasks, borosilicate glass 3.3, Class B. / BRAND
- Volumetric flasks, boro 3.3, class A, with beaded rim / BRAND
- Volumetric flasks, boro 3.3, class A, amber graduations, with PP stopper / BRAND
- Volumetric flasks, boro 3.3, class A, blue graduations, incl. DAkkS calibration certificate / BRAND
- Volumetric Flasks, boro 3.3, class A, blue graduations, with glass stoppers / BRAND
- Volumetric flasks, boro 3.3, class A, blue graduations, with glass stoppers, incl. ISO individual certificate / BRAND
- Volumetric flasks, boro 3.3, class A, blue graduations, with PP stopper, incl. USP batch certificate / BRAND