

Single-channel Pump Heads SP vario for Pump drives Hei-FLOW

- Suitable for most Heidolph tubing due to adjustable roller spacing
- Convex rollers do not squeeze or crush organic cell cultures
- Pump head features ball bearings
- Rollers made of stainless steel and aluminum
- Depending on drive and tubing, the flow rate ranges from 2.0 to 4151 ml per minute

Max. flow rate: 3644 ml/min

Type	Min. flow rate ml/min	PK	Cat. No.
SP vario	0.50	1	9.829 256

Zubehör:
Single-channel Pump Heads SP standard for Pump drives Hei-FLOW / Heidolph Scientific

Alternativen:
Single-channel Pump Heads SP quick for Pump drives Hei-FLOW / Heidolph Scientific