

# Burettes, with straight stopcock, DURAN<sup>&reg;</sup>, class AS, with PTFE spindle

Accuracy class:  
Length scale:  
Stopcock:

AS  
500 mm  
Straight

Nominal capacity ml	Division ml	Tolerance ml	Material stopcock	Colour graduation	With Schellbach stripe	PK	Cat. No.
10	0.02	0.02	PTFE spindle	Blue	Yes	1	9.270 164
25	0.05	0.03	PTFE spindle	Blue	Yes	1	9.270 165
50	0.10	0.05	PTFE spindle	Blue	Yes	1	9.270 166