

Disposable Syringes Omnifix^{®reg}; Solo, 3-piece, Luer slip

- Material: Polypropylene, Plunger backstop: Polyisoprene
- Plunger backstop with smooth double sealing ring for slow aspiration and injection of small sample quantities
- Highly transparent barrel with black graduations (in ml) for ideal readability, permanent marking
- Safe plunger backstop for easy aspiration to the maximum volume
- Inside siliconized
- Luer Slip to put on the needle, centric or eccentric
- Latex- and PVC-free
- Sterile, single packed
- Manufactured according to EN ISO 7886-1

Sterile: Yes
 Package: Single packed

Nominal capacity ml	Material barrel	Material plunger	Colour graduation	Division ml	Nozzle type	Standard	PK	Cat. No.
3	PP	PP	Black	0.1	Luer slip centric	EN ISO 7886-1	100	6.238 513
5	PP	PP	Black	0.2	Luer slip eccentric	EN ISO 7886-1	100	6.301 853
10	PP	PP	Black	0.5	Luer slip eccentric	EN ISO 7886-1	100	7.079 506
20	PP	PP	Black	1	Luer slip eccentric	EN ISO 7886-1	100	6.084 306
30	PP	PP	Black	1	Luer slip eccentric	EN ISO 7886-1	100	6.303 643
50	PP	PP	Black	1	Luer slip eccentric	EN ISO 7886-1	100	4.665 914

Alternativen:

- Fine Dosage Syringes Omnifix[®]-F, 3-piece / B. Braun Deutschland
- Disposable Syringes Omnifix[®] Solo, 3-piece, Luer Lock conus / B. Braun Deutschland