

Drying cabinet OVEN 125 control dry, with glass door

- OVEN control dry with recirculating air up to 300 °C
- Recirculating air drying cabinet with automatic fan switch-off
- Recirculating air adjustable in 10 steps
- Real-time clock
- Temperature selection limited TWB (DIN 12880)
- Illuminated interior
- Cable bushing at the rear for measurings inside the oven
- RS 232 Interface
- Direct saving of data on USB stick
- Interface for additional temperature sensor (display and control possible)
- 2 large TFT Displays

Scope of supply: Oven with 2 shelves and calibration certificate at 160 °C

Temperature:	RT+5 ... 300 °C
Temperature stability:	±0,5 K
Max. no. of inserts:	6
Ext. dim. (WxDxH):	700 x 650 x 825 mm
Int. Dim. (WxDxH):	550 x 450 x 525 mm
Weight:	82 kg
Rated capacity:	2500 W
Plug type:	EU
Power supply:	220 ... 240 V, 50/60 Hz
IP code:	IP 20

Type	Nominal capacity l	Min. temperature °C	Max. temperature °C	With timer	Material interior	Ext. width mm	Ext. dep mm
OVEN 125 control dry	125	+5	300	Yes	Stainless steel	700	650

Zubehör:
Shelves / IKA-Werke