

# High-power homogeniser T 65 digital ULTRA-TURRAX<sup>®</sup>

Extremely powerful dispersing instrument for volumes from 2 to 50 l (H<sub>2</sub>O) for typical pilot plant stations with digital speed display. It offers a wide speed range from 1000 to 9500 rpm that enables users to work at high circumferential speeds even with small rotor diameters. Three rotor-stator configurations guarantee a wide range of applications. Reproducible operations due to constant speed even with changes in viscosities.

- Digital speed display
- Electronic speed control
- Electronic overload protection
- Cables and plugs not included in the delivery
- Plug-in connectors facilitate exchange of dispersing elements
- Stainless steel dispersing elements for use under vacuum and pressure
- Error code display

Working volume:	2 ... 50 l
Speed:	1,000 ... 9,500 1/min
Max. Viscosity:	5000 mPa·s
Noise level:	75 db(A)
Display:	Digital
With timer:	No
Dimensions (WxDxH):	300 x 400 x 420 mm
Weight:	35 kg
Rated capacity:	2600 W
Power supply:	400 V·3/N/PE, 50/60 Hz
IP code:	IP 54

Type	Min. working volume l	Max. working volume l	Min. speed 1/min	Max. speed 1/min	V n 3
T 65 digital ULTRA-TURRAX <sup>®</sup>	2	50	1000	9500	

Zubehör:  
Dispersing tools for Ultra-Turrax<sup>®</sup> T65 basic/digital / IKA-Werke