

Microplates for Cell Culture BRANDplates® cellGrade™

Standard plate for the cultivation of adherent cell cultures. PS-surface with different chemical groups, like e.g. carboxyl and hydroxyl groups, that are freely accessible. Surface is hydrophilic compared with non-treated PS. Incl. one lid per package.

- Sterile
- Manufactured under ISO class 7 cleanroom conditions
- Orange or embossed alphanumeric coding
- Free from endotoxins (< 0.01 EU/ml), DNase, DNA, RNase, non-cytotoxic (according to ISO 10993)
- ANSI/SLAS format

Material: PS
Sterile: Yes
With lid: Yes

No. of wells	Nominal capacity µl	Colour	Bottom shape	RNase/DNase-free	Pyrogen free	Package	PK	Cat. No.
96	330	Transparent	U-shape	Yes	Yes	Single packed	50	4.000 329
96	360	Transparent	V-shape	Yes	Yes	Single packed	50	4.000 330
96	350	Transparent	Flat	Yes	Yes	Single packed	50	4.000 331
96	350	White	Flat	Yes	Yes	Single packed	50	4.000 333
96	350	Black	Flat	Yes	Yes	Single packed	50	4.000 334
96	330	White	Flat	Yes	Yes	Single packed	50	4.000 335
96	330	White	Flat	Yes	Yes	1 x Bag of 5	5	6.272 015
96	330	Black	Flat	Yes	Yes	Single packed	50	4.000 336
96	330	Black	Flat	Yes	Yes	1 x Bag of 5	5	6.272 016
384	100	Transparent	Flat	Yes	Yes	Single packed	50	4.000 337
384	100	White	Flat	Yes	Yes	Single packed	50	4.000 338
384	100	Black	Flat	Yes	Yes	Single packed	50	4.000 339
384	120	White	Flat	Yes	Yes	Single packed	50	4.000 343
384	120	White	Flat	Yes	Yes	1 x Bag of 5	5	6.272 017
384	120	Black	Flat	Yes	Yes	Single packed	50	4.000 344
384	120	Black	Flat	Yes	Yes	1 x Bag of 5	5	6.272 018

Alternativen:

Microplates for Cell Culture BRANDplates® cellGrade™ plus / BRAND
Microplates for Cell Culture BRANDplates® cellGrade™ premium / BRAND
Microplates for Cell Culture BRANDplates® inertGrade™ / BRAND