

Heating Bath HB Digital

Compact and ergonomic heating safety bath for all tempering media.

- Ergonomic handles
- Key button to lock heating bath temperature
- Clear, black-and-white display
- Permanent set-point/actual temperature display
- 3 different safety operating modes
- USB interface

Temperature:	RT ... 180 °C
Dim. Int. (Ø x H):	225 x 130 mm
Ext. dim. (WxDxH):	330 x 325 x 190 mm
Weight:	3.9 kg
Heating capacity:	1350 W
Plug type:	EU
Power supply:	220 ... 240 V, 50/60 Hz
IP code:	IP 21

Type	Nominal capacity l	Max. temperature °C	Display	Operation	Int. Ø mm	Height int. mm	Ext. width mm	Ext. depth mm	Ext. height mm
HB digital	4	180	Digital	Rotary knob	225	130	330	325	190