

Automated cell counter ADAM™ CellT

ADAM™ CellT complies 21 CFR part 11 which a regulation about electronic records and signatures for use in cGMP facilities. When 21 CFR Part 11 mode is activated, data cannot be modified by any user.

- Built-in tablet
- Less than 25 sec to get results
- Automated image analysis
- WiFi (USB Dongle)
- Incl. AccuChip Kit
- 4 TB extra storage for data backup

Dimensions (WxDxH): 276 x 277 x 270 mm
Weight: 7 kg
Power supply: 12 VċDC

Type	Max. volume sample μl	Width mm	Depth mm	Height mm	Type of auxiliary energy	PK	Cat. No.
ADAM™ CellT	13	276	277	270	Mains connection	1	4.684 380

Alternativen:
Automated cell counter ADAM™ MC2 / NanoEn Tek Inc.