

Chamber racks with drawers for ultra-low temperature freezers UF V, stainless steel

Set of 3 chamber racks with drawers. For safe storage of sensitive samples in standard cryoboxes.

Height box:

50 mm

Type of rack:

Drawers

No. of compartments	Arrangement compartments	Width mm	Depth mm	Height mm	Dimensions (WxDxH)	Material	With cryobox
20	5 x 4	140	562	278	140 x 562 x 278 mm	Stainless steel	No
20	5 x 4	140	562	278	140 x 562 x 278 mm	Stainless steel	Yes
24	6 x 4	140	562	334	140 x 562 x 334 mm	Stainless steel	Yes

Alternativen:

Chamber racks with drawers for ultra-low temperature freezers UF V, stainless steel, set of 4 / BINDER