

# Volumetric flasks, boro 3.3, class A, blue graduations, incl. DAkkS calibration certificate

All Blaubrand<sup>®</sup> volumetric flasks are supplied with batch number and one batch certificate per packing unit. They are also available with individual certificates upon request.

- DIN EN ISO 1042
- DE-M marking
- Ring marks and inscriptions in highly contrasting blue enamel
- Calibrated to contain (IN)
- With PP- / PE-stoppers

Standard: DIN EN ISO 1042

Nominal capacity ml	Tolerance ml	Ground size	Accuracy class	Material	Material closure	Type of closure	PK	Cat. No.
25	&plusmn;0.06	NS12/21	A	Boro 3.3	PP	Stopper	1	6.266 579
50	&plusmn;0.06	NS12/21	A	Boro 3.3	PP	Stopper	1	6.266 580
100	&plusmn;0.1	NS12/21	A	Boro 3.3	PP	Stopper	1	6.266 581
1000	&plusmn;0.4	NS24/29	A	Boro 3.3	PP	Stopper	1	6.266 583
2000	&plusmn;0.6	NS29/32	A	Boro 3.3	PP	Stopper	1	6.266 584
10000	&plusmn;2.0	NS45/40	A	Boro 3.3	PE	Stopper	1	6.239 218

Alternativen:  
Volumetric flasks, boro 3.3, class A, blue graduations / BRAND